TRANSPERINEAL NEEDLE GUIDE

REFERENCE GUIDE

For use with
Exact Imaging
EV29L™
transducers

EXACT
IMAGING

NEEDLE GUIDE

This device provides a mechanical means for performing needle / instrument guided procedures with the use of the diagnostic ultrasound endocavity transducer, EV29L.

TRANSPERINEAL NEEDLE GUIDE KIT COMPONENTS

- One (1) EV29L Sterile Transperineal Needle Guide
- One (1) Sterile CIV-FLEX™ Transducer Cover Sheath
- One (1) Sterile Ultrasound Transmission Gel
- Two (2) Orange Kraton Bands

CAUTION


Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING


- Before use, you should be trained in ultrasonography. For complete instructions and warnings related to the use of the needle guide with the EV29L transducer, see your system's user guide.
- See your system's user guide for a glossary of symbols used on the needle guide package.
- Use of a needle guide on an endocavity transducer may increase the risk of vaginal or rectal abrasion.
- For illustration purposes only, transducer may be shown without a transducer cover. Always place a cover over transducer to protect patients and user from cross-contamination.
- Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.
- Disposable components are packaged sterile and are single-use only. Do not use if integrity of packaging is violated or if expiration date has passed.
- Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.
- Do not use if needle guide is damaged or does not fit properly.
- Before use, perform needle path verification to verify ultrasound system and needle guide relationships.
- Use care when inserting needle to avoid bending.

USING TRANSPERINEAL NEEDLE GUIDE


1. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
2. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
3. Secure cover with enclosed bands.
4. Inspect cover to ensure there are no holes or tears.


5. Using proper sterile technique, align locating features of guide body with locating features of transducer.

6. Snap guide body into transducer to secure.


7. Secure latch to bottom of guide body.


SELECTING NEEDLE CHANNEL

1. Once probe is in desired position, slide guide tower into place using grips.

2. Select needle channel and insert needle.

NEEDLE PATH VERIFICATION

- Before use, perform needle path verification to verify ultrasound system and needle guide relationships.
- See transducer user guide for proper procedure.
- Use care when inserting needle to avoid bending.

NOTE: Use appropriate needle length to reach target area.

DISPOSAL

WARNING

- *Dispose of single-use components as infectious waste.*
- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*
- *See your system's user guide for reprocessing transducer between use.*

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Exact Imaging Inc.
Emergo Consulting (UK) Limited c/o Cr360 - UL International