

Material Compatibility List

for Exact Imaging ExactVu™ Transducers

Material Compatibility List (MCL): Listed Chemicals do not degrade transducer performance or safety over 2,500 reprocessing cycles.

NOTE: Microbiological efficacy of the chemicals on the MCL has not been evaluated by Exact Imaging. Please refer to chemical manufacturers' claims.

Type of Chemical	Chemical Name	Manufacturer	Material Compatibility	
			ExactVu™ Endocavity Transducers	
			EV9C	EV29L
Cleaners	Cidezyme	Advanced Sterilization Products	✓	✓
	Enzol	Advanced Sterilization Products	✓	✓
	Aniosyme X3	Laboratories Anios	NC	✓
	Aniosyme DD1	Laboratories Anios	NC	✓
	Revital-Ox Cleaner	Steris	NC	✓
High Level Disinfectants	Aldahol 1.8	Olympus America	NT	✓
	Anioxyde 1000	Laboratories Anios	NC	✓
	Cidex Activated Dialdehyde	Advanced Sterilization Products	✓	✓

	Cidex OPA	Advanced Sterilization Products		
	Metricide Activated Dialdehyde	Metrex Research Inc		
	Revital-Ox Resert	Steris	NC	
	RelyOn PeraSafe	Dupont	NT	

Wipes/Sprays	Acryl-des	Schulke		
	Cleanisept Wipes	Dr. Schumacher GmbH		
	Cleanisept Wipes Forte	Dr. Schumacher GmbH		
	Incidin Foam	Ecolab	NT	
	Tristel Trio	Tristel		
	Tristel Duo	Tristel		
	Surfa'safe	Laboratories Anios	NT	
	Wip'Anios	Laboratories Anios	NT	
	Clinell Universal Wipes	GAMA Healthcare	NT	

Legend:

NC = Not Compatible

NT = Not Tested

Refer to IFU (Care, Cleaning, and Use Guide for EV29L/EV9C) for further instructions

Over time, cleaning and disinfection may cause discoloration.
Discoloration does not affect the performance of the transducer.

Soaking Instructions for ExactVu Transducers:

The images below show the *soak level* for the endocavity transducers.




Figure 1: EV29L Soak Level




Figure 2: EV9C Soak Level

CAUTION

EN-C25



Do not immerse the transducer beyond its soak level.

Rev 1.8

For any further questions about *Cleaning and Disinfection*
with the Exact Imaging transducers,
please contact support@exactimaging.com.